

RODIN[®] BOND

Dental Adhesive System


THE FIRST UNIVERSAL
BONDING AGENT

DESIGNED FOR

3D PRINTED
RESTORATIONS

ALSO DELIVERS

- *High Shear Bond Strength for conventional zirconia and composite restorations*
- *Maximum adhesive strength when coupled with SensiGuard desensitizer, ensuring patient comfort*



*Data from an Independent Research Institute**

RODIN[®] BOND Dental Adhesive System



Rodin Bond demonstrates excellent performance even under high-impact and challenging posterior restorations, with patients who previously experienced post-operative sensitivity reporting no discomfort.

-Dr. Raihan Nazir, DDS, MSc (FFDRCSI)

Product Features

- Unique nanohybrid matrix technology provides exceptional bond strength for 3d printed restorative applications
- Highest adhesion strength for conventional restoratives including Zirconia and composite restorations*
- Efficacy on direct and indirect restorations is further enhanced by SensiGuard™, a dual-function primer and desensitizer
- Flip-top cap design provides “one-drop” flow of adhesive while innovative Unidose system reduces exposure to air and moisture, preserving the bond agent’s effectiveness

*As verified by an independent research institute.




Multiple Delivery Options



Bottle: economic design, precise drop dispensing eliminates waste



Unit Dose: convenience of single-use with precise amount, reduces risk of cross-contamination




RODIN[®] BOND Dental Adhesive System **WORKFLOW**



SENSI GUARD™

When combined with a dental adhesive such as Rodin® Bond, SensiGuard™ acts as both a desensitizer and a primer. Its unique formulation revitalizes collapsed collagen fibers, providing unparalleled monomer penetration, creating exceptional shear bond strength and steadfast adhesion to dentin tubules.

- **Patient Comfort:** Significant reduction of post-operation sensitivity in patients
- **Anti-Bacterial:** Potent antibacterial qualities reduce caries formation beneath restorations
- **Visibility:** Distinct yellow tint provides unparalleled visibility during procedures
- **Universal Compatibility:** Seamless integration with a wide range of dental adhesives and workflows




Longitudinal slices of a sample, subjected to SensiGuard™ treatment, were examined at magnifications of x5,000. The circled region indicates an obstructed area*.

*Data on file.



Scan QR Code to Visit Rodin-3D.com to Learn More about Rodin 3D Resins for Printed Restorations.

RODIN[®] BOND Dental Adhesive System

	Item#	Product Description
	471101	Rodin Bond Dental Adhesive & SensiGuard Desensitizer Intro Kit <i>1 x 5 mL Rodin Bond Bottle</i> <i>1 x 5 mL SensiGuard Bottle</i> <i>2 x 1.2 mL 38% Rodin EtchPro Syringe</i> <i>100 x Micro Brush regular brush</i> <i>20 x Mixing Wells, 4 holes</i> <i>20 x Dispensing tips for PacEtch Pro</i> <i>1 x Technique Guide</i>
	481102	Rodin Bond Dental Adhesive System Unit Dose Kit <i>50 x Rodin Bond 0.1 mL Unit Dose</i> <i>50 x SensiGuard 0.1 mL Unit Dose</i> <i>1 x Technique Guide</i> <i>1 x Instruction Booklet</i>
	471102	Rodin Bond Dental Adhesive Refill, 5mL Bottle
	471103	Rodin Bond Dental Adhesive Unit Dose Refill, 100/pk
	481101	SensiGuard Desensitizer Refill, 5ml Bottle
	481103	SensiGuard Desensitizer Unit Dose Refill, 100/pk
	135	Rodin EtchPro 38% Phosphoric Acid Etch <i>4 x 1.2 mL syringes, 8 x Dispensing Tips</i>
	135E	Econo pack: 36 x 1.2 mL syringes
	135C	Clinic pack: 144 x 1.2 mL syringes
	231	25 Ga tips, Etch, Blue, 100/pk
	901	Disposable MicroPro Applicators Standard 2.0mm, 400/pk
	902	Disposable MicroPro Applicators Fine 1.5mm, 400/pk
	903	Disposable MicroPro Applicators Ultra fine 1.0mm, 400/pk
	904	Disposable MicroPro Applicators Tapered, 400/pk
	906	Mixing Well, 4-Well, 200/pk