

RODIN[®] InteCem

Intermediate Restorative Material

- ✓ Provides sedative like properties for hypersensitive tooth pulp
- ✓ High strength and abrasion resistance with low solubility
- ✓ Acts as a proper thermal insulator
- ✓ Ideal as a base under resin-free restorative materials
- ✓ Remains in place up to one year



pac-dent.com




RODIN[®] InteCem

Intermediate Restorative Material


Rodin InteCem Intermediate Restorative Material is designed for intermediate restorations that can remain in place for up to one year. It features a unique composition of polymer-reinforced zinc oxide-eugenol that leverages eugenol's sedative properties to relieve hypersensitive tooth pulp and provide effective thermal insulation for enhanced patient comfort.

Rodin InteCem is an ideal base for resin-free cements and restorative materials, including amalgam, inlays, and onlays. It boasts high strength comparable to zinc phosphate cements, with added benefits of high abrasion resistance, optimal sealing properties, and low solubility.


Net Setting Time (min)



Working Time (min)



Comprehensive Strength (MPa)



BUY 3, GET 1 FREE!

452201	Rodin InteCem Kit: 1 x 38 g Powder & 1 x 14 mL Liquid	NET: \$45.135 Reg. \$60.18
---------------	---	--------------------------------------



pac-dent.com



© 2024 Pac-Dent, Inc. All rights reserved.

Pac-Dent[®]
Passion for Excellence