

PacEndo™
VITALICE™
360 GREEN
 Pulp Vitality Refrigerant Spray

COLDER CLEANER GREENER



PacEndo™ Vital-Ice™ Green 360 Pulp Vitality Spray is an ultra-low-GWP cold spray for endodontic thermal testing, providing a fast, consistent cold stimulus without refrigeration. Delivers an ultra-cold temperature of -26.2° C – much colder than ice or ethyl chloride for proper testing.

Vital-Ice™ 360 Green is formulated with a next-generation propellant featuring a short atmospheric lifetime, offering clinicians a more environmentally friendly choice that supports reliable chairside responses. With a Global Warming Potential <1 and zero chlorofluorocarbons (CFCs), Vital-Ice™ 360 Green is the ideal choice for responsible dental professionals.


Compare to CanalPro® Endo-Ice® Green Cold Spray

- Environmentally friendly low GWP <1, contains no chlorofluorocarbons (CFCs)
- Extremely cold temperature at -26.2° C
- Produces a rapid test response
- No refrigeration needed
- Spray nozzle allows for precise application

CanalPro® and Endo-Ice® are registered trademarks of Coltène/Whaledent

Net: 29.99 MSRP: 39.99

911101 | PacEndo™ Vital-Ice™ 360 Green, 1 x 6oz Bottle



VITAL ICE 360 GREEN	VITAL ICE
✓ GWP <1	✗ GWP <1
✓ No CFCs	✗ No CFCs
✓ No Refrigeration Needed	✓ No Refrigeration Needed
✓ Precision Spray Nozzle	✓ Precision Spray Nozzle

BUY 3 AND GET 1 FREE

To qualify and receive your promotional goods please submit a valid proof of purchase to freegoods@pac-dent.com. Invoices must be submitted via email to qualify. Cannot Be Combined With Any Other Offer. No Cash Value. Same or Lesser Item(s) Applied. Offers are valid only at participating dealers in the continental US, excluding AK and HI. Promotions are subject to change or cancellation without notice. Please allow 6-8 weeks for delivery.