


PAC-DENT PROPHY ANGLE FAMILY

Disposable Propphy Angle Comparison Chart



Prophy Angle	AntiSplatr® Flex Slim	AntiSplatr®	Oscillating AntiSplatr Mini	Mini 90®	Mighty Mini 90®	Mini Ergo®	ProAngle® EZ	ProAngle® Ergo	FUNimal™	ProAngle® Plus	ProAngle® Plus with Brush Cup	ProAngle®
SKU	212101 212102 212103 212104	LSFAS-100 LSFAF-100 LSFAS-500 LSFAF-500	ASMS-100 ASMF-100 ASMB-100	\$163-90-144 \$164-90-144 \$163-90-500 \$164-90-500	M162-90-144 M163-90-144 M163-90-500	\$163-144 \$164-144 \$163-500 \$164-500	EZMS-100 EZES-500 EZSS-100 EZMF-100 EZEF-500 EZSF-100 EZMS-500 EZES-100 EZSS-500 EZMF-500 EZEF-100 EZSF-500	168-144 168-500 168S-144 168S-500 188-144 188-500 198T-144 198T-500	166-100A 167-100P 166-500A 167-100R	249-144 249-50 249S-144 249S-500 250-144 250-500 251T-144 251T-500	238-100 239-100	247-144 248-144 247-500 248-500
MSRP	\$0.83 - 0.85/Each	\$0.78 - 0.79/Each	\$1.05 - 1.16/Each	\$0.89 - 0.93/Each	\$0.77 - 0.78/Each	\$0.89 - 0.93/Each	\$0.76 - 0.80/Each	\$0.66/Each	\$0.68 - 1.05/Each	\$0.45 - 0.50/Each	\$0.84/Each	\$0.42 - 0.43/Each
Body Size	Slim	Regular	Mini	Mini	Mini	Mini	Extra Long	Regular	Regular	Regular	Regular	Regular
Body Type	90 Degree / Contra Angle	90 Degree	90 Degree	90 Degree	90 Degree	Contra Angle	90 Degree	Contra Angle	90 Degree	90 Degree	90 Degree	90 Degree
Description	◆ Ideal contour adaptation and interproximal cleaning ◆ Flex cup minimizes cup slipping during procedures while reducing paste splatter ◆ 42% smaller body and neck design for enhanced access and visibility ◆ Splatter fender minimizes the accumulation of paste and bodily fluids ◆ Ergo configuration reduces hand fatigue	◆ Splatter fender removes accumulated paste and bodily fluids ◆ Helps protect against the risk of cross-contamination ◆ Spiral blades allow for optimum paste retention ◆ Optimal cup lip and flare allows for greater subgingival access ◆ Reinforces existing infection control protocols	◆ Oscillating spin reduces splatter ◆ Lowered heat generation allows for increased contact ◆ Ultra-compact design expands visibility and reach ◆ Unique cup design provides optimal flare around tooth contours ◆ Quieter and smoother than traditional unidirectional DPAs	◆ Compact head and slimmer profile allow superior visibility and handling ◆ Improved flare for better subgingival access ◆ Patented beveled gear design eliminates vibration and ensures smooth, quiet operation ◆ Pre-sterilized with gamma radiation to reduce the risk of cross-contamination	◆ Compact, powerful head for increased visibility and improved access ◆ Improved flare for better subgingival access ◆ Patented beveled gear design delivers quiet and smooth power for enhanced patient comfort	◆ Ergonomically designed with the optimal amount of head tilt to reduce hand and wrist fatigue ◆ 20% shorter head and 25% slimmer neck profile for increased visibility and improved access ◆ Patented beveled gear design eliminates vibration and ensures smooth, quiet operation ◆ Pre-sterilized with gamma radiation to reduce the risk of cross-contamination	◆ Reduce the weight of the handpiece, which reduces hand and wrist fatigue for dentists and hygienists ◆ Eliminates the costs of cleaning and disinfection for the nose cone ◆ Patented beveled gear design, the angle runs smoothly with no overheating	◆ Ergonomic design reduces hand and arm micro-traumas ◆ The gearing structure of ProAngle Ergo makes it much quieter than other contra-angle DPAs ◆ Torque Cup minimizes splatter, improves efficiency by 15%, and offers better interproximal cleaning	◆ Designed with a compact head, shorter neck, and a slimmer profile that makes polishing procedures for young patients with smaller teeth easier and more engaging than ever ◆ These FUN angles will reduce your young patients' anxiety while offering optimal performance	◆ Reinforced beveled gears to ensure only one angle is required for each procedure. ◆ Eliminates vibration for less noise, chatter, and heat ◆ Guaranteed no cup dropping	◆ Beveled gears are precisely engineered so that their pitch surfaces are mounted securely and firmly when engaged	◆ Patented prophy angle that employs the same design as ProAngle Plus to ensure each prophy procedure only requires one angle to finish ◆ Utilizes the same beveled gear design as ProAngle Plus for assured quality and performance
Indications	Patients with small mouths, deep occlusal anatomy, or orthodontic appliances requiring precision cleaning	General Prophy, high-risk aerosol infection control cases, pediatric or geriatric patient who require extra cleanliness	Pediatric, geriatric, or special needs patients requiring gentle and low-vibration, low heat polishing	Children, patients with limited mouth opening, and detailed polishing around orthodontic brackets	Patients with heavy stain or calculus, as well as areas requiring precision polishing	High-volume hygiene practices and clinicians with repetitive strain concerns	General prophylaxis procedures, especially in busy dental practices focusing on efficiency	Hygienists performing multiple cleanings daily, patients with tight interproximal spaces	Pediatric dentistry and nervous or anxious children who benefits from a fun, less intimidating experience	Standard hygiene procedures, heavy workload settings where reliability and comfort are key	Heavy stain removal, orthodontic bracket cleaning, and prophylaxis in patients with deep fissures	Routine dental cleaning, general prophylaxis, and patients with moderate plaque buildup
Configuration	200/pk, 1200/pk	100/pk, 500/pk	100/pk	144/pk, 500/pk	144/pk, 500/pk	144/pk, 500/pk	100/pk, 500/pk	144/pk, 500/pk	100/pk, 500/pk	144/pk, 500/pk	100/pk	144/pk, 500/pk
Compatible Handpieces	Doriot	Doriot	Doriot	Doriot	Doriot	Doriot	E-Type, M-Type, S-Type	Doriot	Doriot	Doriot	Doriot	Doriot
Cup Options	Flex	Soft, Firm	Soft, Firm, Tapered Brush	Soft, Firm	Long Soft, Soft	Soft, Firm	Soft, Firm	Soft, Super Soft, Firm, Torque	Pedo Soft, Soft	Soft, Super Soft, Firm, Torque	Tapered Brush, Flat Brush	Soft, Firm
Ergonomically design	Yes	No	No	No	No	Yes	No	Yes	No	No	No	No
Optimal Flare	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	N/A	Yes
Sterile	No	No	No	Yes	No	Yes	No	No	No	No	No	No
Latex-Free	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes



Scan QR code to request free samples