

EzAim X-Ray Solution

Product Catalog




TABLE OF CONTENTS

EzAim® Family	1
Color Guide	2
EzAim® Fit	3
EzAim® Grip	4
EzAim® Snap	4
EzAim® Defend 3-in-1	5
EzAim® ComforTab	6
EzAim®	7
X-Ray Accessories	8
Product Sku Guild	9


EzAim® Family

X-Ray Positioning System & Accessories




EzAim® Fit

Universal Adjustable Holder
Page: 3




EzAim® Grip

Ratcheting Sensor Holder
Page: 4




EzAim® Snap

Sensor Holder
Page: 4




EzAim® ComforTab

Adhesive Sensor Holder
Page: 6



EzAim® Defend 3-in-1

Collimator Positioning System
Page: 5



EzAim®

X-Ray Positioning System
Page: 7


Find Your Perfect Fit

Our diverse lineup is engineered to meet every practice's specific requirements:


- **Reusable Excellence:** The Fit, Grip, and Snap products offer autoclavable durability and feature a cushioned bite tab to significantly reduce patient discomfort.
- **Disposable Convenience:** For high-volume efficiency, the ComforTab and Individually Cut Holders provide a hygienic, single-use solution with strong adhesive for reliable positioning.
- **Universal Compatibility:** Most EzAim holders are designed to fit all sensor types on the market and work seamlessly with your existing aiming rings and positioning arms.

EzAim® Color Guide


Anterior




Bitewing



Posterior



Endo



Sensor Holder Comparasion




Product Names	Fit	Grip	Snap	ComforTab	Individually Cut Holder	Original EzAim Holder
Cushioned Bite Tab	✓	✓	✓	✓	✗	✗
Individually Cut & Packed	✓	✓	✓	✓	✓	✗
Strong Adhesive	✗	✓	✗	✓	✓	✓
Non-Slip Surface	✓	✓	✓	✓	✗	✗
Reduces Patient Discomfort	✓	✓	✓	✓	✗	✗
Universal	✓	✓	✗	✓	✓	✓
Reusable	✓	✗	✓	✗	✗	✗
Item #	131135, 131136, 131137, 131138, 131139, 131140	131148, 131149	131150, 131151, 131152, 131153	131115, 131121, 131122, 131123, 131124, 131125	131116, 131117, 131118, 131119, 131120	131109, 131110, 131111, 131112, 131113

EzAim® Fit


Universal Adjustable Holder



Universal Adjustable Holder



Universal Fit for Every Sensor



"Added Cushion for Comfort in Every Bite"

Compare to Rinn® XCP-DS Fit® Sensor Holders

The new EzAim Fit Universal Sensor Holders are an advanced solution that simplifies the imaging process with fewer parts to keep track of. The color-coded expandable sensor holders fit most x-ray sensors, sizes 1, 1.5, and 2 on the market. Features a unique, cushioned silicone bite area and cord clip for added patient comfort. Compatible with existing aiming rings and positioning arms.

- Expandable universal design holder fits all digital sensors
- Features cushioned silicone bite tab for extra patient comfort
- Includes built in cord clip for seamless cord management
- Autoclavable design for low cost per use

BUY 3, GET 1 of the Same FREE

EzAim® Grip



Universal Ratcheting Holder



131149

Compare to Rinn® Uni-Grip® 360 Sensor Holders

The EzAim Grip Disposable Universal Sensor Holder is engineered for full-mouth efficiency, featuring a precision ratcheting tab that rotates to capture posterior, anterior, and bitewing positions — all with a single holder. Designed to integrate seamlessly with the EzAim Defend 3-in-1 System, it eliminates the need for multiple holders during a full mouth series, saving time and reducing cross-contamination risk. A cushioned silicone bite tab improves patient comfort, and perforated tabs allow for quick, mess-free removal after use.

- Disposable ratcheting tab rotates for full-mouth imaging with one holder
- Cushioned silicone bite tab enhances patient comfort across all positions
- Accommodates sensor sizes #1, #1.5, and #2 for universal compatibility
- Compact, smooth-edged tab design minimizes patient discomfort & oral tissue contact

BUY 3, GET 1 of the Same FREE


EzAim® Snap



Autoclavable Sensor Holder



131151



131152



131153



131150

Compare to Instrumentarium KaVo IXS™ Digital Sensor Holders


The New EzAim Snap Autoclavable X-Ray Sensors are designed to be used with the KaVo IXS™ and Gendex™ GXS-700 digital sensor holders. Features a cushioned silicone bite tab for added patient comfort and built in cord clip for seamless cable management.

- Features cushioned silicone bitetab for extra patient comfort
- Color-coded holders for size 1 and 2 sensors
- Autoclavable design for low cost per use

BUY 3, GET 1 of the Same FREE

EzAim® Defend 3-in-1

X-Ray Collimator Positioning System



3-in-1 Simplicity for Smarter Imaging

Compare to Rinn® XCP-ORA®

The EzAim Defend 3-in-1 Collimator Positioning System features a compact and lightweight design that eliminates extra parts needed for dental imaging, streamlining the imaging process and protecting patients from unnecessary radiation exposure.

The easily adjustable, color-coded positioning ring allows for quick and accurate alignment, making it ideal for bitewing, anterior, and posterior exposures. Compatible with all EzAim bite blocks, the system works seamlessly with digital sensors, phosphor plates, and traditional film, ensuring adaptability to various clinical setups. Fully autoclavable arm and ring.


- Built-in collimator reduces radiation & enhances diagnostics
- 3-in-1 color-coded ring for bitewing, anterior & posterior shots
- Lightweight, compact design—no extra parts needed
- Compatible with all EzAim bite blocks for digital sensors, phosphor plates, & film
- Autoclavable arm and ring



Red Bitewing





BUY 3, GET 1 of the Same FREE




EzAim® ComforTab™

Adhesive Cushioned Sensor Holders



Universal Adhesive Fit for Every Sensor



Individually Cut!

Compare to Schick AimRight™ and Xios Plus™


The EzAim ComforTabs are the newest addition to the EzAim X-Ray Positioning System, designed to elevate patient comfort without compromising imaging accuracy. Featuring a soft, cushioned bite tab, ComforTabs offer a more comfortable experience during radiographic procedures while creating a non-slip surface that holds sensors securely in place. The gentle yet stable design reduces patient discomfort and minimizes sensor movement for more consistent, high-quality images.

- **Cushioned Tab for Extra Patient Comfort**
- **Individually Cut and Packaged**
- **Strong Adhesive Allows for Reliable Positioning**


BUY 3, GET 1 of the Same FREE

EzAim®

X-Ray Positioning System



Try out the NEW Short Arm Option



Sensor Holders Compared to Xios Plus

Arms and Rings Compared to Dentsply Rinn® XCP & XCP Short Arm

The EzAim X-Ray Positioning System features disposable, self-adhesive sensor holders with peel-and-stick adhesive technology that firmly secures digital sensors to holders. The holders are color coded for easy identification and fits all types of sensors on the market. EzAim holders can be used with existing aiming rings and positioning arms. Patient-friendly design eliminates extra dimensions for setup. Can be managed with hemostats. Available with matching, autoclavable positioning arms and aiming rings.

- **Single-use adhesive holder that optimizes hygienic workflow**
- **Color coded for easy identification**
- **Saves time and eliminates reprocessing**
- **Available for anterior, posterior, endo tab, and bite tab**

BUY 3, GET 1 of the Same FREE

X-Ray Accessories

Armor™ X-Ray Protection



Armor™ Digital X-Ray Sensor Sheaths

These specially designed sheaths with paper-backing are custom fitted for each digital X-ray sensor type to ensure maximum patient comfort and ease of use. The sheaths' films are made of specially formulated high quality blend of resins that are durable and scratch-resistant, protecting your investment, while saving money on labor.

Each package contains **500 sleeves**.

BUY 3, GET 1 of the Same FREE



Armor™ Phosphor Storage Plate Barrier Envelopes

Using Y-notch design that ensures the envelope tears open in a single direction to reduce contamination. Armor Barrier Envelopes provide an extra layer of defense, allowing Phosphor Storage Plate (PSP) sensors to be easily reused. The black opaque back of the envelope prevents ambient light exposure, preventing image distortion.

Compare to Air Techniques' ScanX® Barrier Envelopes

BUY 3, GET 1 of the Same FREE

SoftEdge™ Cushions X-Ray Film and Sensor Cushions



SoftEdge Cushions features disposable, peel-and-stick adhesive cushions that easily adhere to digital sensors, phosphor plates, and film packets. SoftEdge Cushions create a soft, comfortable edge to protect and safeguard soft tissue from abrasions while allowing patients to bite harder on the positioning system. Provides technicians with improved diagnostic results due to proper periapical positioning, actively reducing the need for retakes.

Compare to Crosstex® Edge-Ease® Cushions

SKU: 131301 | Blue Cushions, Large
131302 | White Cushions, Small

BUY 3, GET 1 of the Same FREE

Bite Blocks



SureHold™ Disposable X-Ray Film Holder

SureHold™ Disposable X-Ray Film Holders effectively and comfortably holds size #0, #1, or #2 film for anterior and posterior exposures for XCP® and BAI film positioning systems. Convenient scored front easily snaps off for pedo applications and imaging.

SKU: 131202



PropBlock™ Disposable Bite Blocks

PropBlock™ Disposable Bite effectively and comfortably holds size #0, #1, or #2 film for use with the XCP® film positioning system for vertical or horizontal exposures. Convenient scored front easily snaps off for pedo applications and imaging.

SKU: 131203



EasyProp™ Disposable Mouth Props

EasyProp™ Disposable Mouth Props helps patients hold the desired angle, without worrying about movement or angle disruption. The integrated hole design allows technicians to thread a strand of dental floss when needed.

SKU: 131201

BUY 3, GET 1 of the Same FREE

Try out the NEW EzAim Products

BUY 3, GET 1 of the Same FREE!

ITEM#	DESCRIPTION	COMPARED TO	QUANTITY
	Comfortab		
131121	Blue Anterior ComforTab	Xios Plus 6176510	50/pack
131122	White Anterior ComforTab	Schick AimRight B1073054	50/pack
131123	Blue Bitewing ComforTab	Schick AimRight B1073052	50/pack
131124	Red Bitewing ComforTab	Xios Plus 6176536	50/pack
131125	Yellow Posterior ComforTab	Schick AimRight B1073055	50/pack
131115	Red SoftBite Tab Sensor Holder	Schick AimRight B1073051	100/pack
	Defend 3-in-1 Positioning System		
131127	EzAim Defend 3-in-1 Positioning System, 1/set	XCP-ORA 550771	1/set
131128	EzAim Defend 3-in-1 Positioning System, 2/set	XCP-ORA 550772	2/set
131129	EzAim Defend 3-in-1 Positioning Ring	XCP-ORA 550773	1/pack
	Positioning Arms		
131102	Red Bitewing Positioning Arm	XCP/BAI 540927	1/pack
131126	Blue Anterior/Yellow Posterior Positioning Arm	XCP-ORA 550774	1/pack
131131	Blue Short Anterior Positioning Arm	XCP Rin Short Arm 540871	1/pack
131132	Red Short Bitewing Positioning Arm	XCP Rin Short Arm 540873	1/pack
131133	Yellow Short Posterior Positioning Arm	XCP Rin Short Arm 540872	1/pack
131134	Green Short Endo Positioning Arm	XCP Rin Short Arm 540874	1/pack
131130	EzAim Complete Kit: ComforTab with Defend 3-in-1		
	1x Defend 3-in-1 Positioning Ring and 2x Positioning Arms, 60 ComforTab Holders (10 of each: Blue Anterior, White Anterior, Red Bitewing, Blue Bitewing, Yellow Posterior, and Bitetab)		
	EzAim Fit NEW		
131135	Anterior Thin Biteblock	XCP DS FIT 559906	2/pack
131136	Anterior Biteblock	XCP DS FIT 559901	2/pack
131137	Horizontal Bitewing Biteblock	XCP DS FIT 559903	2/pack
131138	Vertical Bitewing Biteblock	XCP DS FIT 559904	2/pack
131139	Posterior Biteblock	XCP DS FIT 559902	2/pack
131140	Endodontic Biteblock	XCP DS FIT 559905	2/pack
131141	Periapical Size 0	XCP DS FIT 559907	2/pack
131142	Fit Complete Kit: 1x Defend Ring & arm set, 2x each biteblock, 1x endo ring, 1x endo arm, 50x universal sleeves		
131143	Fit Hygiene Kit:1x Defend Ring & arm set, 2x each anterior, posterior, VBW, HBW biteblocks, 50x universal sleeves		
131144	Fit Elastic short band refill (posterior, VBW, Endo)	XCP DS FIT 559913	6/pack
131145	Fit Elastic long band refill (anterior, HBW)	XCP DS FIT 559914	6/pack
	EzAim Grip NEW		
131148	EzAim Grip Universal Sensor Holder Intro Kit, 1x Defend ring & arm, with 100 universal sensor sleeves	XCP Uni-Grip 550252	100/pack
131149	EzAim Grip Universal Sensor holders refill kit, with 50 universal sensor sleeves	XCP Uni-Grip 550052	50/pack
	EzAim Snap NEW		
131150	Snap Size 1 Anterior Holder - For IXS Sensors	Dexis IXS 1.013.2905	1/pack
131151	Snap Size 1 HBW Holder - For IXS Sensors	Dexis IXS 1.013.2906	1/pack
131152	Snap Size 1 VBW Holder - For IXS Sensors	Dexis IXS 1.013.2907	1/pack
131153	Snap Size 1 Posterior Holder - For IXS Sensors	Dexis IXS 1.013.2904	1/pack
131154	Snap Size 2 Anterior Holder - For IXS Sensors	Dexis IXS 1.013.2910	1/pack
131155	Snap Size 2 HBW Holder - For IXS Sensors	Dexis IXS 1.013.2911	1/pack
131156	Snap Size 2 VBW Holder - For IXS Sensors	Dexis IXS 1.013.2912	1/pack
131157	Snap Size 2 Posterior Holder - For IXS Sensors	Dexis IXS 1.013.2909	1/pack
131173	EzAim Snap Size 1 Intro Kit: 1x Defend Set, 2x HBW, VBW, Posterior, Anterior holders		
131174	EzAim Snap Size 2 Intro Kit: 1x Defend Set, 2x HBW, VBW, Posterior, Anterior holders		



*Trademark of Dentsply Sirona



670 Endeavor Circle Brea, CA 92821 | (909) 839-0888 |
info@pac-dent.com | www.pac-dent.com